

Anti-HPV33 Major capsid protein L1 Recombinant Antibody

Catalog No: #V1099

Package Size: #V1099-1 50ug #V1099-2 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Anti-HPV33 Major capsid protein L1 Recombinant Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	VS2206
Isotype	IgG1
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA-cap, Neutralization
Species Reactivity	Human papillomavirus 33 (HPV33)
Conjugates	Unconjugated
Target Name	L1, Major capsid protein L1
Accession No.	P06416
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>95% as determined by SDS-PAGE.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.