

Anti-HPV35 E7/Protein E7 Recombinant Antibody

Catalog No: #V1101



Package Size: #V1101-1 50ug #V1101-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-HPV35 E7/Protein E7 Recombinant Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	VS1007
Isotype	IgG2a, kappa
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, IHC
Species Reactivity	Human papillomavirus type 35 (HPV35)
Conjugates	Unconjugated
Accession No.	P27230
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.

Application Details

Purity:>95% purity as determined by SDS-PAGE.ELISA: 1:2000-1:20000

Images



SDS-PAGE for Anti-HPV35 E7/Protein E7 Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.