

Anti-HPV45 Major capsid protein L1 Recombinant Antibody

Catalog No: #V1107

Package Size: #V1107-1 50ug #V1107-2 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Anti-HPV45 Major capsid protein L1 Recombinant Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	VS2172
Isotype	IgG
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA
Species Reactivity	Human papillomavirus 45
Conjugates	Unconjugated
Target Name	L1, Major capsid protein L1
Accession No.	P36741
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>95% as determined by SDS-PAGE.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.