

Anti-HPV45 Major capsid protein L1 Recombinant Antibody

Catalog No: #V1108

Package Size: #V1108-1 50ug #V1108-2 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Anti-HPV45 Major capsid protein L1 Recombinant Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	VS2090
Isotype	IgG
Purification	Protein A/G purified from cell culture supernatant.
Species Reactivity	Human papillomavirus 45
Conjugates	Unconjugated
Target Name	L1, Major capsid protein L1
Accession No.	P36741
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>95% as determined by SDS-PAGE.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.