

Anti-HPV52 Major capsid protein L1 Recombinant Antibody

Catalog No: #V1110

Package Size: #V1110-1 50ug #V1110-2 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Anti-HPV52 Major capsid protein L1 Recombinant Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	VS1427
Isotype	IgG2b
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA-cap, Neutralization
Species Reactivity	Human papillomavirus 52 (HPV52)
Conjugates	Unconjugated
Target Name	L1, Major capsid protein L1
Accession No.	Q05138
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>95% as determined by SDS-PAGE.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.