

Anti-PEG Recombinant Antibody

Catalog No: #V2498



Package Size: #V2498-1 50ug #V2498-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-PEG Recombinant Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	EG1003
Isotype	IgG1, kappa
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, FCM, WB
Species Reactivity	General
Conjugates	Unconjugated
Target Name	PEG, Poly(ethylene glycol), Polyethylene glycol, 6.3, 6-3
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Application Details

Purity:>95% as determined by SDS-PAGE.ELISA: 1:2000-1:20000; WB: 1:500-1:2000

Images



SDS-PAGE for Anti-PEG Antibody .

Note: This product is for in vitro research use only and is not intended for use in humans or animals.