

Anti-PEG Recombinant Antibody

Catalog No: #V2503



Package Size: #V2503-1 50ug #V2503-2 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Anti-PEG Recombinant Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	EG1014
Isotype	IgG
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, WB
Species Reactivity	General
Conjugates	Unconjugated
Target Name	PEG, Poly(ethylene glycol), Polyethylene glycol
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>95% as determined by SDS-PAGE.ELISA: 1:2000-1:20000; WB: 1:500-1:2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.