

Anti-ssDNA/dsDNA/Z-DNA Recombinant Antibody

Catalog No: #V2644



Package Size: #V2644-1 50ug #V2644-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-ssDNA/dsDNA/Z-DNA Recombinant Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	NA1178
Isotype	IgG2a, kappa
Purification	Protein A/G purified from cell culture supernatant.
Applications	DB, ELISA
Conjugates	Unconjugated
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Application Details

Purity:>95% as determined by SDS-PAGE.ELISA: 1:2000-1:20000

Images



SDS-PAGE for Anti-ssDNA/dsDNA/Z-DNA Antibody .

Note: This product is for in vitro research use only and is not intended for use in humans or animals.