

Anti-dsDNA Recombinant Antibody

Catalog No: #V2677



Package Size: #V2677-1 50ug #V2677-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-dsDNA Recombinant Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	NA1201
Isotype	IgG2a, kappa
Purification	Protein A/G purified from cell culture supernatant.
Applications	DB, ELISA, IF
Conjugates	Unconjugated
Formulation	0.94 mMB Citric acid monohydrate, 16.7mMB Trisodium citrate dihydrate, 6.2%B sucrose, pHB 6.5, 0.01%B Tween 80
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Application Details

Purity:>95% as determined by SDS-PAGE.ELISA: 1:2000-1:20000

Images



SDS-PAGE for Anti-dsDNA Antibody .

Note: This product is for in vitro research use only and is not intended for use in humans or animals.