

Anti-DNA/RNA G-quadruplex Structures Recombinant Antibody

Catalog No: #V2706



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Package Size: #V2706-1 50ug #V2706-2 100ug


Description

Product Name	Anti-DNA/RNA G-quadruplex Structures Recombinant Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	NA1053
Isotype	IgG2a, kappa
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, EMSA, ICC, IF, ChIP (ChIP-seq), WB
Conjugates	Unconjugated
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Application Details

Purity:>95% as determined by SDS-PAGE.ELISA: 1:2000-1:20000; WB: 1:500-1:2000

Images



SDS-PAGE for Anti-DNA/RNA G-quadruplex Structures Antibody .

Note: This product is for in vitro research use only and is not intended for use in humans or animals.