

# Anti-Triplex DNA Recombinant Antibody

Catalog No: #V2707



Package Size: #V2707-1 50ug #V2707-2 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|              |                                                                                                                                                                 |
|--------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name | Anti-Triplex DNA Recombinant Antibody                                                                                                                           |
| Host Species | Mouse                                                                                                                                                           |
| Clonality    | Monoclonal                                                                                                                                                      |
| Clone No.    | NA1078                                                                                                                                                          |
| Isotype      | IgG2a, kappa                                                                                                                                                    |
| Purification | Protein A/G purified from cell culture supernatant.                                                                                                             |
| Applications | DB, ELISA, IF                                                                                                                                                   |
| Conjugates   | Unconjugated                                                                                                                                                    |
| Formulation  | 0.01M PBS pH 7.4                                                                                                                                                |
| Storage      | Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks).<br>Store at -20°C 12 months. Store at -80°C long term. |

## Application Details

Purity:>95% as determined by SDS-PAGE.ELISA: 1:2000-1:20000

## Images



SDS-PAGE for Anti-Triplex DNA Antibody .

Note: This product is for in vitro research use only and is not intended for use in humans or animals.