

Anti-DNA G-quadruplex Structures Recombinant Nanobody

Catalog No: #V2710

Package Size: #V2710-1 50ug #V2710-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-DNA G-quadruplex Structures Recombinant Nanobody
Host Species	Alpaca
Clonality	Monoclonal
Clone No.	NA1029
Isotype	VHH-His
Purification	Purified by Nickel column.
Applications	ELISA, EMSA, IF
Conjugates	Unconjugated
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>95% as determined by SDS-PAGE.ELISA: 1:2000-1:20000

Images



SDS-PAGE for Anti-DNA G-quadruplex Structures Nanobody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.