

Recombinant RNA Crystalization Chaperone Fab

Catalog No: #V2731

Package Size: #V2731-1 50ug #V2731-2 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant RNA Crystalization Chaperone Fab
Host Species	Human
Clonality	Monoclonal
Clone No.	NA1079
Isotype	Fab - G1 - kappa
Purification	Protein A/G purified from cell culture supernatant.
Applications	Crystalization
Conjugates	Unconjugated
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>90% as determined by SDS-PAGE.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.