

Anti-REGN5381 Reference Antibody

Catalog No: #V2732




Package Size: #V2732-1 50ug #V2732-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Anti-REGN5381 Reference Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	AD2066
Isotype	IgG
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, Bioactivity:B FACS, FunctionalB assay, ResearchB inB vivo
Conjugates	Unconjugated
Target Name	Anti-REGN5381 Antibody
Formulation	0.01M PBS pH 7.4

Images



SDS-PAGE for Research Grade Anti-REGN5381 Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.