

Anti-FMC63 scFv Monoclonal Recombinant Antibody

Catalog No: #V2733




Package Size: #V2733-1 50ug #V2733-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-FMC63 scFv Monoclonal Recombinant Antibody
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	AD2075
Isotype	IgG
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, FCM
Conjugates	Unconjugated
Target Name	FMC63
Formulation	0.01M PBS pH 7.4

Images



Flow-cytometry using Anti-FMC63 scFv Monoclonal Antibody. FMC63 Transfected CHO cells were stained with an irrelevant antibody (Blue Histogram) or an Anti-FMC63 scFv monoclonal antibody. (Green Histogram) at a concentration of 5 µg/ml for 30 mins at RT. After washing, bound antibody was detected using a FITC conjugated goat anti-rabbit antibody and cells analysed on a NovoCyte Flow Cytometer.



SDS-PAGE for Anti-FMC63 scFv Monoclonal Antibody .

Note: This product is for in vitro research use only and is not intended for use in humans or animals.