

Anti-Trastuzumab Idiotypic Recombinant Antibody

Catalog No: #V2736




Package Size: #V2736-1 50ug #V2736-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Anti-Trastuzumab Idiotypic Recombinant Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	AD2011
Isotype	IgG2a, kappa
Purification	Protein A/G purified from cell culture supernatant.
Applications	Blocking, ELISA
Conjugates	Unconjugated
Formulation	0.01M PBS pH 7.4

Images



SDS-PAGE for Anti-Trastuzumab Idiotypic Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.