

Anti-FMC63 scFv Monoclonal Recombinant Antibody

Catalog No: #V2737




Package Size: #V2737-1 50ug #V2737-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Anti-FMC63 scFv Monoclonal Recombinant Antibody
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	AD2200
Isotype	IgG
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA
Conjugates	Unconjugated
Target Name	FMC63
Formulation	0.01M PBS pH 7.4

Images



SDS-PAGE for Anti-FMC63 scFv Monoclonal Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.