

Anti-Abatacept Recombinant Antibody

Catalog No: #V2749




Package Size: #V2749-1 50ug #V2749-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-Abatacept Recombinant Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	AD2023
Isotype	IgG1, kappa
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA
Conjugates	Unconjugated
Target Name	Cytotoxic T-lymphocyte protein 4, Cytotoxic T-lymphocyte-associated antigen 4, CTLA-4, CD152, CTLA4, CD152
Formulation	0.01M PBS pH 7.4

Images



Detects Abatacept in indirect ELISA



SDS-PAGE for Anti-Abatacept Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.