

Anti-Etanercept Non-Neutralizing Recombinant Antibody

Catalog No: #V2750


Package Size: #V2750-1 50ug #V2750-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Anti-Etanercept Non-Neutralizing Recombinant Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	AD2126
Isotype	IgG2a, kappa
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA
Conjugates	Unconjugated
Target Name	Etanercept
Formulation	0.01M PBS pH 7.4

Images



SDS-PAGE for Anti-Etanercept Non-Neutralizing Recombinant Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.