

Anti-Adalimumab Non-Neutralizing Recombinant Antibody

Catalog No: #V2751

Package Size: #V2751-1 50ug #V2751-2 100ug




Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Anti-Adalimumab Non-Neutralizing Recombinant Antibody
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	AD2106
Isotype	IgG
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA
Conjugates	Unconjugated
Formulation	0.01M PBS pH 7.4

Images



Detects Adalimumab in indirect ELISA

Note: This product is for in vitro research use only and is not intended for use in humans or animals.