

Anti-Adalimumab Neutralizing Recombinant Antibody

Catalog No: #V2753



Package Size: #V2753-1 50ug #V2753-2 100ug


Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Anti-Adalimumab Neutralizing Recombinant Antibody
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	AD2028
Isotype	IgG
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA
Conjugates	Unconjugated
Formulation	0.01M PBS pH 7.4


Images

Anti-Adalimumab Neutralizing Antibody (AD2028)
Blocking with TNFa



Blocking Adalimumab binds with TNFa

Anti-Adalimumab Neutralizing Antibody (AD2028)
binds with Adalimumab



Detects Adalimumab in indirect ELISA

Note: This product is for in vitro research use only and is not intended for use in humans or animals.