

Anti-Adalimumab Neutralizing Recombinant Antibody

Catalog No: #V2754



Package Size: #V2754-1 50ug #V2754-2 100ug


Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Anti-Adalimumab Neutralizing Recombinant Antibody
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	AD2067
Isotype	IgG
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA
Conjugates	Unconjugated
Formulation	0.01M PBS pH 7.4


Images

Anti-Adalimumab Neutralizing Antibody (AD2067)
binds with Adalimumab



Detects Adalimumab in indirect ELISA

Anti-Adalimumab Neutralizing Antibody (AD2067)
Blocking with TNF α

Blocking Adalimumab binds with TNF α

Note: This product is for in vitro research use only and is not intended for use in humans or animals.