

# Anti-Pembrolizumab Neutralizing Recombinant Antibody

Catalog No: #V2773

Package Size: #V2773-1 50ug #V2773-2 100ug


Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                      |
|-----------------------|------------------------------------------------------|
| Product Name          | Anti-Pembrolizumab Neutralizing Recombinant Antibody |
| Host Species          | Mouse                                                |
| Clonality             | Monoclonal                                           |
| Clone No.             | AD2261                                               |
| Isotype               | IgG1, kappa                                          |
| Purification          | Protein A/G purified from cell culture supernatant.  |
| Applications          | ELISA                                                |
| Immunogen Description | Pembrolizumab                                        |
| Conjugates            | Unconjugated                                         |
| Target Name           | Pembrolizumab, Keytruda, MK-3475, lambrolizumab      |
| Formulation           | 0.01M PBS pH 7.4                                     |


## Images

Anti-Pembrolizumab Neutralizing Antibody (AD2261)  
Blocking with PD1




Blocks the binding of Pembrolizumab to PD1.

Anti-Pembrolizumab Neutralizing Antibody (AD2261)  
binds with Pembrolizumab



Detects Pembrolizumab in indirect ELISA.



Note: This product is for in vitro research use only and is not intended for use in humans or animals.