

Anti-Evolocumab Neutralizing Recombinant Antibody

Catalog No: #V2781




Package Size: #V2781-1 50ug #V2781-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-Evolocumab Neutralizing Recombinant Antibody
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	AD2268
Isotype	IgG
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA
Conjugates	Unconjugated
Formulation	0.01M PBS pH 7.4

Images



Detects Evolocumab in indirect ELISA.



Blocks the binding of PCSK9 to Evolocumab.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.