

Anti-Tocilizumab Neutralizing Recombinant Antibody

Catalog No: #V2792



Package Size: #V2792-1 50ug #V2792-2 100ug


Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Anti-Tocilizumab Neutralizing Recombinant Antibody
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	AD2220
Isotype	IgG
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA
Conjugates	Unconjugated
Formulation	0.01M PBS pH 7.4


Images

Anti-Tocilizumab Neutralizing Antibody (AD2220)
binds with Tocilizumab



Detects Tocilizumab in indirect ELISA.

Anti-Tocilizumab Neutralizing Antibody (AD2220)
Blocking with IL6R



SAA2627 blocks the binding of Tocilizumab to IL-6R.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.