

Anti-Guselkumab Neutralizing Recombinant Antibody

Catalog No: #V2804



Package Size: #V2804-1 50ug #V2804-2 100ug


Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Anti-Guselkumab Neutralizing Recombinant Antibody
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	AD2180
Isotype	IgG
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA
Conjugates	Unconjugated
Formulation	0.01M PBS pH 7.4


Images

Anti-Guselkumab Neutralizing Antibody (AD2180)
Blocking with IL23A



Blocks the binding of Guselkumab to IL23A.

Anti-Guselkumab Neutralizing Antibody (AD2180)
binds with Guselkumab



Detects Guselkumab in indirect ELISA.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.