

# Anti-Guselkumab Neutralizing Recombinant Antibody

Catalog No: #V2805




Package Size: #V2805-1 50ug #V2805-2 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|              |                                                     |
|--------------|-----------------------------------------------------|
| Product Name | Anti-Guselkumab Neutralizing Recombinant Antibody   |
| Host Species | Rabbit                                              |
| Clonality    | Monoclonal                                          |
| Clone No.    | AD2231                                              |
| Isotype      | IgG                                                 |
| Purification | Protein A/G purified from cell culture supernatant. |
| Applications | ELISA                                               |
| Conjugates   | Unconjugated                                        |
| Formulation  | 0.01M PBS pH 7.4                                    |

## Images



Blocks the binding of Guselkumab to IL23A.



Detects Guselkumab in indirect ELISA.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.