

Anti-Ara h 6 Human IgE Recombinant Antibody

Catalog No: #V3029



Package Size: #V3029-1 50ug #V3029-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Anti-Ara h 6 Human IgE Recombinant Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	AG3133
Isotype	IgE
Applications	ELISA
Species Reactivity	Arachis hypogaea (Peanut)
Conjugates	Unconjugated
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Application Details

ELISA: 1:2000-1:20000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.