

Anti-Aspergillus fumigatus mitF/aspF1 Recombinant Antibody

Catalog No: #V3031



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Package Size: #V3031-1 50ug #V3031-2 100ug


Description

Product Name	Anti-Aspergillus fumigatus mitF/aspF1 Recombinant Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	AG3178
Isotype	IgG1
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, WB
Species Reactivity	Neosartorya fumigata (strain ATCC MYA-4609 / Af293 / CBS 101355 / FGSCA1100) (Aspergillus fumigatus)
Conjugates	Unconjugated
Accession No.	P67875
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Application Details

ELISA: 1:2000-1:20000; WB: 1:500-1:2000

Images



SDS-PAGE for Anti-Aspergillus fumigatus mitF/aspF1 Antibody .

Note: This product is for in vitro research use only and is not intended for use in humans or animals.