

Anti-Bet v 1/Allergen Bet v I-A Recombinant Antibody

Catalog No: #V3034



Package Size: #V3034-1 50ug #V3034-2 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	Anti-Bet v 1/Allergen Bet v I-A Recombinant Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	AG3060
Isotype	IgG1
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA
Species Reactivity	Betula pendula (European white birch)
Conjugates	Unconjugated
Accession No.	P15494
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Application Details

ELISA: 1:2000-1:20000

Images

Anti-Bet v 1/Allergen Bet v I-A Antibody (AG3060)
binds with Bet v 1/Allergen Bet v I-A



Detects Recombinant European white birch Bet v 1/Allergen
Bet v I-A Protein, N-His in indirect ELISA.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.