

Anti-Bet v 1/Allergen Bet v I-A Recombinant Antibody

Catalog No: #V3044



Package Size: #V3044-1 50ug #V3044-2 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


Description

|                    |                                                                                                                                                                 |
|--------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name       | Anti-Bet v 1/Allergen Bet v I-A Recombinant Antibody                                                                                                            |
| Host Species       | Human                                                                                                                                                           |
| Clonality          | Monoclonal                                                                                                                                                      |
| Clone No.          | AG3017                                                                                                                                                          |
| Isotype            | VHH-hFc                                                                                                                                                         |
| Purification       | Protein A/G purified from cell culture supernatant.                                                                                                             |
| Applications       | ELISA, WB, SPR                                                                                                                                                  |
| Species Reactivity | Betula pendula (European white birch)                                                                                                                           |
| Conjugates         | Unconjugated                                                                                                                                                    |
| Accession No.      | P15494                                                                                                                                                          |
| Formulation        | 0.01M PBS pH 7.4                                                                                                                                                |
| Storage            | Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks).<br>Store at -20°C 12 months. Store at -80°C long term. |


Application Details

ELISA: 1:2000-1:20000

Images



SDS-PAGE for Anti-Bet v 1/Allergen Bet v I-A Antibody .



Lane 1: recombinant European white birch Bet v 1/Allergen Bet v I-A Lane 2: 1\*loading buffer

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.