

Anti-European white birch Bet v 1/Allergen Bet v I-A Recombinant Nanobody

Catalog No: #V3048

Package Size: #V3048-1 50ug #V3048-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-European white birch Bet v 1/Allergen Bet v I-A Recombinant Nanobody
Host Species	Alpaca
Clonality	Monoclonal
Clone No.	AG3122
Isotype	VHH-8His-Cys-tag
Purification	Purified by Nickel column.
Applications	ELISA, SPR
Species Reactivity	Betula pendula (European white birch)
Conjugates	Unconjugated
Accession No.	P15494
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

ELISA: 1:2000-1:20000

Images



SDS-PAGE for Anti-European white birch Bet v 1/Allergen Bet v I-A Nanobody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.