

Anti-Bet v 1/Allergen Bet v I-A Recombinant Antibody

Catalog No: #V3049

Package Size: #V3049-1 50ug #V3049-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-Bet v 1/Allergen Bet v I-A Recombinant Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	AG3037
Isotype	IgG1, kappa
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, SPR
Species Reactivity	Betula pendula (European white birch)
Conjugates	Unconjugated
Accession No.	P15494
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.


Application Details


ELISA: 1:2000-1:20000

Images




Detects recombinant European white birch Bet v 1/Allergen Bet v I-A in indirect ELISA





SDS-PAGE for Anti-Bet v 1/Allergen Bet v I-A Antibody .



Western blot analysis was performed using anti-Bet v I monoclonal antibody at 1ug/mL on various samples. Lane 1: recombinant European white birch Bet v 1/Allergen Bet v I-A
Lane 2: negative control

Note: This product is for in vitro research use only and is not intended for use in humans or animals.