

Anti-BETVII/Allergen Bet v II Monoclonal Antibody

Catalog No: #V3050




Package Size: #V3050-1 50ug #V3050-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Anti-BETVII/Allergen Bet v II Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	AG3176
Isotype	IgG2a
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, WB
Species Reactivity	Betula pendula (European white birch) (Betula verrucosa)
Immunogen Description	Recombinant BETVII/Allergen Bet v II protein.
Conjugates	Unconjugated
Formulation	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Images



SDS-PAGE for Anti-BETVII/Allergen Bet v II Monoclonal Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.