

Anti-BETVII/Allergen Bet v II Recombinant Antibody

Catalog No: #V3051



Package Size: #V3051-1 50ug #V3051-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-BETVII/Allergen Bet v II Recombinant Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	AG3045
Isotype	IgG1
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, WB
Species Reactivity	Betula pendula (European white birch) (Betula verrucosa)
Conjugates	Unconjugated
Accession No.	P25816
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Application Details

ELISA: 1:2000-1:20000; WB: 1:500-1:2000

Images



SDS-PAGE for Anti-BETVII/Allergen Bet v II Antibody .

Note: This product is for in vitro research use only and is not intended for use in humans or animals.