

Anti-Bovine Serum Albumin/BSA Monoclonal Antibody

Catalog No: #V3057




Package Size: #V3057-1 50ug #V3057-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Anti-Bovine Serum Albumin/BSA Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	AG3187
Isotype	IgG1
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA
Species Reactivity	Bovine
Immunogen Description	Natural BSA.
Conjugates	Unconjugated
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Images



SDS-PAGE for Anti-Bovine Serum Albumin/BSA Monoclonal Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.