

Anti-Bovine Serum Albumin/BSA Recombinant Antibody

Catalog No: #V3060

Package Size: #V3060-1 50ug #V3060-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-Bovine Serum Albumin/BSA Recombinant Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	AG3089
Isotype	IgG1, lambda
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA
Species Reactivity	Bovine
Conjugates	Unconjugated
Accession No.	P02769
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Application Details

ELISA: 1:2000-1:20000

Images



SDS-PAGE for Anti-Bovine Serum Albumin/BSA Antibody .

Note: This product is for in vitro research use only and is not intended for use in humans or animals.