

Anti-Bovine Serum Albumin/BSA Recombinant Antibody , HRP

Catalog No: #V3067



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Package Size: #V3067-1 50ug #V3067-2 100ug


Description

Product Name	Anti-Bovine Serum Albumin/BSA Recombinant Antibody , HRP
Host Species	Human
Clonality	Monoclonal
Clone No.	AG3104
Isotype	IgG1, kappa
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA (Det), WB
Species Reactivity	Bovine
Conjugates	Unconjugated
Accession No.	P02769
Formulation	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Application Details

ELISA: 1:2000-1:20000; WB: 1:500-1:2000

Images



SDS-PAGE for Anti-Bovine Serum Albumin/BSA Antibody , HRP

Note: This product is for in vitro research use only and is not intended for use in humans or animals.