

Anti-Bovine Serum Albumin/BSA Recombinant Antibody , Biotin

Catalog No: #V3068

Package Size: #V3068-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Anti-Bovine Serum Albumin/BSA Recombinant Antibody , Biotin
Host Species	Human
Clonality	Monoclonal
Clone No.	AG3107
Isotype	IgG1, kappa
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA (Det), WB
Species Reactivity	Bovine
Conjugates	Unconjugated
Accession No.	P02769
Formulation	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

ELISA: 1:2000-1:20000; WB: 1:500-1:2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.