

Anti-Can f 1 Human IgE Recombinant Antibody

Catalog No: #V3071

Package Size: #V3071-1 50ug #V3071-2 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Anti-Can f 1 Human IgE Recombinant Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	AG3091
Isotype	IgE
Applications	ELISA
Species Reactivity	Canis lupus familiaris (Dog) (Canis familiaris)
Conjugates	Unconjugated
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Application Details

ELISA: 1:2000-1:20000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.