

Anti-Can f 1 Recombinant Antibody

Catalog No: #V3072

Package Size: #V3072-1 50ug #V3072-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-Can f 1 Recombinant Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	AG3139
Isotype	IgG
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA
Species Reactivity	Canis lupus familiaris (Dog) (Canis familiaris)
Conjugates	Unconjugated
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

ELISA: 1:2000-1:20000

Images



SDS-PAGE for Anti-Can f 1 Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.