

Anti-Dermatophagoides farinae DERF1/DERP1 Recombinant Antibody

Catalog No: #V3079

Package Size: #V3079-1 50ug #V3079-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description


Product Name	Anti-Dermatophagoides farinae DERF1/DERP1 Recombinant Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	AG3130
Isotype	IgG1, kappa
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, WB
Species Reactivity	Dermatophagoides farinae (American house dust mite), Dermatophagoides pteronyssinus (European house dust mite)
Conjugates	Unconjugated
Accession No.	P16311, P08176
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Application Details


ELISA: 1:2000-1:20000; WB: 1:500-1:2000

Images


Anti-Dermatophagoides farinae DERF1/DERP1 Antibody (AG3130) binds with DERF1



Detects Dermatophagoides farinae (American house dust mite), Dermatophagoides pteronyssinus (European house dust mite) DERF1/DERP1 in indirect ELISAs.



SDS-PAGE for Anti-Dermatophagoides farinae
DERF1/DERP1 Antibody .



Western blot analysis was performed using anti-American
house dust mite DERF1 monoclonal antibody at 1ug/mL on
various samples. Lane 1: recombinant American house dust
mite Der f 1/DERF1 Lane 2: negative control

Note: This product is for in vitro research use only and is not intended for use in humans or animals.