

Anti-Dermatophagoides pteronyssinus DERP2/DERPII Recombinant Antibody

Catalog No: #V3099



Package Size: #V3099-1 50ug #V3099-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-Dermatophagoides pteronyssinus DERP2/DERPII Recombinant Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	AG3099
Isotype	IgG1, kappa
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA
Species Reactivity	Dermatophagoides pteronyssinus (European house dust mite)
Conjugates	Unconjugated
Accession No.	P49278
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Application Details

ELISA: 1:2000-1:20000

Images



SDS-PAGE for Anti-Dermatophagoides pteronyssinus DERP2/DERPII Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.