

Anti-Feline Allergen Fel d 1 Reference Antibody

Catalog No: #V3103


Package Size: #V3103-1 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-Feline Allergen Fel d 1 Reference Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	AG3018
Isotype	IgG4
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, Bioactivity:B FACS, FunctionalB assay, ResearchB inB vivo
Species Reactivity	Feline
Conjugates	Unconjugated
Accession No.	P30438, P30440
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Images



Detects *Felis catus* (Cat) (*Felis silvestris catus*) Fel d 1-B/CH2/Allergen FdI in indirect ELISAs.



SDS-PAGE for Research Grade Anti-Cat Allergen Fel d 1 Antibody .

Note: This product is for in vitro research use only and is not intended for use in humans or animals.