

# Anti-CH1&CH2/Allergen Fel d I-A&B Recombinant Antibody

Catalog No: #V3104

Package Size: #V3104-1 50ug #V3104-2 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                    |                                                                                                                                                                 |
|--------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name       | Anti-CH1&CH2/Allergen Fel d I-A&B Recombinant Antibody                                                                                                          |
| Host Species       | Human                                                                                                                                                           |
| Clonality          | Monoclonal                                                                                                                                                      |
| Clone No.          | AG3040                                                                                                                                                          |
| Isotype            | IgG4                                                                                                                                                            |
| Purification       | Protein A/G purified from cell culture supernatant.                                                                                                             |
| Applications       | ELISA, SPR                                                                                                                                                      |
| Species Reactivity | Cat                                                                                                                                                             |
| Conjugates         | Unconjugated                                                                                                                                                    |
| Accession No.      | P30438, P30440                                                                                                                                                  |
| Formulation        | 0.01M PBS pH 7.4                                                                                                                                                |
| Storage            | Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks).<br>Store at -20°C 12 months. Store at -80°C long term. |

## Application Details

ELISA: 1:2000-1:20000

## Images



SDS-PAGE for Anti-CH1&CH2/Allergen Fel d I-A&B Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.