

Anti-Denosumab Polyclonal Antibody

Catalog No: #V6849



Package Size: #V6849-1 50ug #V6849-2 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Anti-Denosumab Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by antigen affinity column.
Applications	ELISA
Immunogen Description	Denosumab
Conjugates	Unconjugated
Target Name	Rabbit polyclonal to Prolia.
Formulation	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.