

Anti-Sarilumab Polyclonal Antibody

Catalog No: #V6858

Package Size: #V6858-1 50ug #V6858-2 100ug



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Anti-Sarilumab Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by antigen affinity column.
Applications	ELISA
Immunogen Description	Sarilumab
Conjugates	Unconjugated
Target Name	Rabbit polyclonal to Kevzara.
Formulation	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.