

Anti-Tocilizumab Polyclonal Antibody

Catalog No: #V6859



Package Size: #V6859-1 50ug #V6859-2 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Anti-Tocilizumab Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by antigen affinity column.
Applications	ELISA
Immunogen Description	Tocilizumab
Conjugates	Unconjugated
Target Name	Rabbit polyclonal to Actemra.
Formulation	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.