

Anti-Rituximab Polyclonal Antibody

Catalog No: #V6861

Package Size: #V6861-1 50ug #V6861-2 100ug



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Anti-Rituximab Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by antigen affinity column.
Applications	ELISA
Immunogen Description	Rituximab
Conjugates	Unconjugated
Target Name	Rabbit polyclonal to Rituxan.
Formulation	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.