

Anti-Pembrolizumab Polyclonal Antibody

Catalog No: #V6868

Package Size: #V6868-1 50ug #V6868-2 100ug



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Anti-Pembrolizumab Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by antigen affinity column.
Applications	ELISA
Immunogen Description	Pembrolizumab
Conjugates	Unconjugated
Target Name	Rabbit polyclonal to Keytruda.
Formulation	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.